



Office of the Medical Director
Noel Wagner, MD, NREMT-P
1000 Houghton Ave
Saginaw, MI 48602
(989) 583-7940
Fax (989) 583-7941
www.SaginawValleyEMS.org

June 26, 2014

MEMORANDUM #2014-03

FROM: Office of the Medical Director
Saginaw Valley Medical Control Authority (SVMCA)

SUBJECT: PHARMACY RECERTIFICATION – June 27, 2014

The Greater Bay Area EMS Medication System will be recertifying all EMS drug boxes and kits on **June 27, 2014**. All boxes and bags will need to be turned into a participating pharmacy for recertification during the morning hours beginning around 07:30 am.

Recertified boxes/kits will have an **orange** tag which will signify that they are good through October, 2014

CHANGES:

Naloxone

The amount of naloxone in the drug boxes is being doubled. As you know, there has been an increase in the number and severity of narcotic overdoses. Since the overdosed drugs have a half-life greater than naloxone; it may be necessary to give higher doses to maintain adequate respirations in some patients.

Glucagon will be moved into the same spot as the Solu-Medrol and the additional naloxone will take its place.

Nitroglycerin Drips

The nitroglycerin drips will remain removed from the boxes for at least the next four months. Although supply levels are returning to normal, potential shortages could reoccur.

Controlled Substances Kits

We have purchased new Pelican kits to replace the canvas kits. The new kits will be sealed through a hole that has been drilled into the kit. The number of the box is engraved on the right-upper corner of the kit. We will be replacing the old canvas kits periodically until all of the old kits have been removed.

There are no changes to either the ALS or BLS soft kits.

As always, if you have any question please feel free to contact our office.

A handwritten signature in black ink, appearing to read 'Eric Snidersich'.

Eric Snidersich, BS EMT-P
EMS Manager, SVMCA

Attachment: Inventory List (June 2014)



| MEDICATION BOX INVENTORY - JUNE 2014 | | | |
|--------------------------------------|-------------------------------------|--|--|
| TRAY 1 | | TRAY 3 | |
| 2 | Clave Adapter | | |
| 2 | Atomizer | | |
| 3 | IV Start Kit | | |
| 2 | Saline Flush | 10 ml | pre-filled syringe |
| 3 | Safety IV Catheter | 14 GA x 1.75" | Insyte Autoguard |
| 3 | Safety IV Catheter | 16 GA x 1.16" | Insyte Autoguard |
| 3 | Safety IV Catheter | 18 GA x 1.16" | Insyte Autoguard |
| 3 | Safety IV Catheter | 20 GA x 1.16" | Insyte Autoguard |
| 3 | Safety IV Catheter | 22 GA x 1" | Insyte Autoguard |
| 3 | Safety IV Catheter | 24 GA x 1" | Insyte Autoguard |
| TRAY 2 | | MEDICATIONS - BOTTOM OF BOX | |
| 2 | Diphenhydramine | 50 mg/1 ml | vial |
| 1 | Methylprednisolone | 125 mg/2 ml | vial |
| 1 | Glucagon injection | 1 mg | diluent/powder |
| 4 | Naloxone | 2 mg/2 ml | pre-filled syringe |
| 2 | Acetaminophen (Tylenol) | 160 mg/5 ml | unit-dose cup |
| 1 | Oral Syringe | 5 ml | |
| 2 | Ondansetron (Zofran) | 4 mg/2 ml | vial |
| 3 | Safety Needle | 19 GA x 1 1/2" | |
| 3 | Safety Needle | 21 GA x 1 1/2" | |
| 3 | Safety Needle | 23 GA x 5/8" | |
| 2 | Ipratropium (Atrovent) | 0.5 mg/2.5 ml | unit-dose vial |
| 2 | Albuterol (Proventil) | 2.5 mg/3 ml | unit-dose vial |
| 2 | Magnesium Sulfate | 1 Gram/2 ml | vial |
| 1 | Epinephrine 1:1000 | 30 mg /30 ml | vial |
| 1 | Nitroglycerin | 0.4 mg tab (#25) | bottle |
| 4 | Aspirin | 81 mg chew tab | unit-dose packet |
| 1 | Red plastic reseal lock | | |
| 2 | Vasopressin | 20 u/1 ml | vial |
| 3 | Amiodarone | 150 mg/3 ml | vial |
| 1 | Dopamine | 400 mg/5 ml | vial |
| 5 | Adenosine | 6 mg/2 ml | vial |
| 1 | Bacteriostatic NAACL | 30 ml | vial |
| 3 | Atropine Sulfate | 1 mg/10 ml | pre-filled syringe |
| 2 | Lidocaine | 100 mg/5ml | pre-filled syringe |
| 1 | Calcium Chloride 10% | 1 Gram/10 ml | pre-filled syringe |
| 2 | Syringe | 1 ml | |
| 2 | Syringe | 3 ml | Luer Lock |
| 2 | Syringe | 10 ml | Luer Lock |
| 1 | Syringe | 20 ml | Luer Lock |
| 1 | Syringe | 30 or 35 ml | Luer Lock |
| TRAY 2 | | MEDICATIONS - BOTTOM OF BOX | |
| 1 | Dextrose 5% in Water | 25 Grams/50 ml | pre-filled syringe |
| 8 | Epinephrine 1:10,000 | 1 mg/10 ml | pre-filled syringe |
| 2 | Sodium Bicarbonate 8.4% | 50 MEq/50 ml | pre-filled syringe |
| TRAY 2 | | IV SOLUTIONS - BOTTOM OF BOX | |
| 1 | Dextrose 5% in Water | 250 ml | |
| 2 | Normal Saline | 500 ml | |
| 1 | Dextrose 5% in Water 50ml | 50 ml | 1 bag minimum. May be more. |
| TRAY 2 | | IV SETS - BOTTOM OF BOX | |
| 3 | Primary Macro Drip Tubing | | Clave port |
| 2 | Primary Micro Drip Tubing | | Clave port |
| 1 | Primary Pump Tubing | | Plum pump |
| TRAY 2 | | MISCELLANEOUS - BOTTOM OF BOX | |
| 1 | Nitroglycerin (temporarily removed) | 50 mg/250 ml | pre-mixed container |
| 1 | Paper Bag | | |
| 1 | Nebulizer | | |
| 1 | 4-way Stopcock Extension | | |
| 1 | Illinois Illiac Needle | 16 ga. (8881245164) | Kendall-Monoject |
| 1 | Illinois Illiac Needle | 15 ga. (DIN1515X) | Baxter |
| TRAY 2 | | RED BAG-CONTROLLED SUBSTANCES (SEALED WITH RED) | |
| 1 | Carpuject Injector | | |
| 2 | Diazepam | 10 mg/2 ml | carpuject or vial |
| 2 | Fentanyl | 100 mcg/2 ml | carpuject or vial |
| 2 | Morphine | 10 mg/1 ml | carpuject, vial, or pre-filled syringe |
| 2 | Midazolam | 10 mg/2 ml | vial |

| ALS MINI-KIT - JUNE 2014 | | | |
|--------------------------|---------------------------|------------------|--------------------|
| 2 | Clave Adapter | | |
| 1 | IV Start Kit | | |
| 2 | Saline Flush Syr | 10 ml | |
| 2 | Safety IV Catheter | 18 GA x 1.16" | Insyte Autoguard |
| 2 | Safety IV Catheter | 20 GA x 1.16" | Insyte Autoguard |
| 2 | Safety IV Catheter | 22 GA x 1" | Insyte Autoguard |
| 1 | Syringe | 3 ml | Luer Lock |
| 1 | Safety Needle | 19 GA x 1 1/2" | |
| 1 | Ondansetron (Zofran) | 4 mg/2 ml | vial |
| 2 | Ipratropium (Atrovent) | 0.5 mg/2.5 ml | unit-dose vial |
| 2 | Albuterol (Proventil) | 2.5 mg/3 ml | unit-dose vial |
| 1 | Nitroglycerin | 0.4 mg tab (#25) | bottle |
| 4 | Aspirin | 81 mg chew tab | unit-dose packet |
| 1 | 50% Dextrose | 25 Grams/50 ml | pre-filled syringe |
| 1 | Normal Saline | 250 ml | |
| 1 | Primary Macro Drip Tubing | | Clave port |
| 1 | Red plastic reseal lock | | |
| 1 | Paper Bag | | |

| BLS KIT - JUNE 2014 | | | |
|---------------------|-----------------------|----------------|------------------|
| 1 | Epi-Pen | 0.3 mg | Auto Injector |
| 1 | Epi-Pen Jr. | 0.15 mg | Auto Injector |
| 1 | Nebulizer | | |
| 4 | Aspirin | 81 mg chew tab | unit-dose packet |
| 2 | Albuterol (Proventil) | 2.5 mg/3 ml | unit-dose vial |